

SUBMISSION FEES & INSTRUCTIONS

CANINE			FELINE		
SERVICE REQUESTED	SPECIMEN	FEE	SERVICE REQUESTED	SPECIMEN	FEE
Anaplasma phagocytophilum titer (IFA)	serum	\$28.00	FAVN Rabies titer (for export)	Serum	\$75.00
Autogenous Wart Vaccine	warts in formalin	\$90.00	Feline calicivirus antibody titer (IFA)	serum	\$15.00
Canine Snap 4DX Plus	serum	\$22.00	Feline coronavirus (FCoV) antibody titer	serum	\$20.00
CAV-2 antibody - IgG only (vaccine titer)	serum	\$15.00	Feline panleukopenia- IgG only (vaccine titer)	serum	\$15.00
CDV antibody - IgG only (vaccine titer)	serum	\$15.00	Feline panleukopenia titer (IgG & IgM)	serum	\$24.00
CDV antibody IgG & IgM (IFA)	serum / CSF	\$24.00	FeLV and FIV Combo	EDTA blood, plasma, serum	\$28.00
CHV antibody titer (IFA)	serum	\$15.00	FeLV IFA antigen stain	EDTA blood, bone marrow	\$20.00
CHV antigen stain (FA)	vaginal swabs	\$18.00	FHV antibody titer (IFA)	serum	\$15.00
CPV antibody - IgG only (vaccine titer)	serum	\$15.00	FHV antigen stain (FA)	nasal/ conjunctival slides	\$18.00
CPV antibody IgG & IgM (IFA)	serum	\$24.00	Rabies antibody titer (RFFIT)	serum	\$50.00
CPV antigen - feces	feces	\$22.00	Toxoplasma antibody titer IgG & IgM (IFA)	serum / CSF	\$25.00
Ehrlichia canis titer (IFA)	serum	\$24.00	Virus isolation (Call lab before submitting)		\$48.00
FAVN Rabies titer (for export)	Serum	\$75.00	BOVINE		
Lyme titer (IFA)	serum	\$23.00	SERVICE REQUESTED	SPECIMEN	FEE
Rabies antibody titer (RFFIT)	serum	\$50.00	Autogenous Wart Vaccine	warts in formalin	\$90.00
RMSF titer (IFA)	serum	\$24.00	EQUINE		
Toxoplasma antibody titer IgG & IgM (IFA)	serum / CSF	\$25.00	SERVICE REQUESTED	SPECIMEN	FEE
Virus isolation (Call lab before submitting)		\$48.00	Autogenous Wart Vaccine	warts in formalin	\$90.00
			Lyme titer (IFA)	serum	\$23.00
			Rabies antibody titer (RFFIT)	serum	\$50.00

SAMPLE COLLECTION

SERUM: When serum is required, DO NOT send whole blood. DO NOT send extremely hemolyzed and/or lipemic samples.

SLIDES: DO NOT fix slides before sending, air dry only.

FIXED TISSUE: Please place wart tissue in neutral buffered formalin. The volume of formalin used should be AT LEAST 5X the weight of the tissue, i.e. ¼ cup of wart tissue should be in 1 ¼ cup of buffered formalin.

SWABS: Place swabs in sterile PBS or viral transport media. DO NOT use any of the bacteriological gel transport systems.

SHIPPING INSTRUCTIONS

Serum, CSF, swabs, EDTA blood, and effusions should be packed in a leak proof container with absorbent material. This package should be placed inside a second container on ice (gel packs). Samples should be clearly labeled. Overnight shipping is recommended. Ship via UPS or Federal Express. Please do not ship samples on Friday.

Slides and fixed tissues can be sent without ice packs.

CREDIT CARD PAYMENTS ARE AVAILABLE ONLINE: <https://www.aub.ie/payinvoice>